

What Bottles are Compatible with Hygeia Flanges?

Picture	Manufacturer and Bottle Type	Connects without Adapter	Connects with Hygeia Adapter	Not Compatible
	Adiri NxGen Nurser			X
	Ameda Breastmilk Storage Bottles, 4oz		X	
	Avent - Natural Feeding, 4 oz			X
	Philips Avent - Glass			X
	Philips Avent - Classic, 9 oz			X
	Babies R Us Purely Simple Feeding Bottles, 5oz			X
	Born Free Classic - Plastic, 5oz			X
	Breast Flow - The First Years			X
	Comotomo Silicone Bottle			X
	Dr. Brown's Glass Bottle, 4 oz		X	
	Dr. Brown's Standard Natural Flow		X	
	Dr. Brown's Wide Neck			X
	EvenFlo Breastmilk storage containers, 5oz		X	
	EvenFlo BPA Free Classic Glass/Plastic, 4oz		X	
	Gerber First Essentials Bottles, 5oz		X	
	Lansinoh mOmma Feeding Bottle, 250 mL			X
	Lansinoh mOmma Bottle with Natural Wave, 5oz	X		
	LifeFactory Glass Bottles		X	
	Medela - 80mL Storage Container		X	

Picture	Manufacturer and Bottle Type	Connects without Adapter	Connects with Hygeia Adapter	Not Compatible
	Medela - 150mL Storage Container		X	
	Medela - Breast Milk Bottle		X	
	Medela - Calma 8oz Set		X	
	Munchkin Mighty Grip			X
	Nuby - Standard Neck		X	
	Nuby - Wide Neck			X
	Nuk Orthodontic Wide mouth			X
	Playtex - Drop In Nurser			X
	Playtex -Vent Aire BPA Free Wide Bottle			X
	Similac Bottle		X	
	Simplisse Breastmilk Collection Bottles, 4oz		X	
	ThinkBaby - 5oz and 9oz	X		
	Tommee Tippee - Closer to Nature			X